

## Kinase Assay

<u>Lysis buffer</u> (final conc.):	for 20 ml:
• 20 mM Tris/HCl pH7.5	400 $\mu$ l 1 M
• 150 mM NaCl	600 $\mu$ l 5 M
• 25 mM $\beta$ -glycerophosphate	500 $\mu$ l 1 M
• 2 mM EDTA	80 $\mu$ l 0.5 M
• 2 mM pyrophosphate	400 $\mu$ l 0.1 M
• 1 mM orthovanadate	200 $\mu$ l 0.1 M
• 1% Triton X-100	2 ml 10%
• 1 mM DTT	20 $\mu$ l 1 M
• 1 mM NaF	20 $\mu$ l 1 M
• A. dest.	15.8 ml

Protease Inhibitors: added before use (Leupeptin, Pepstatin, Pefa-Block) according to stock

<u>Kinase buffer</u> (final conc.):	for 20 ml:
• 20 mM Tris/HCl pH7.5	400 $\mu$ l 1 M
• 20 mM $\beta$ -glycerophosphate	400 $\mu$ l 1 M
• 100 $\mu$ M orthovanadate	20 $\mu$ l 0.1 M
• 10 mM MgCl <sub>2</sub>	200 $\mu$ l 1 M
• 50 mM NaCl	200 $\mu$ l 5 M
• 1 mM DTT	20 $\mu$ l 1 M
• 50 $\mu$ M ATP	50 $\mu$ l 20 mM
• 1 mM NaF	20 $\mu$ l 1 M
• A. dest.	18.7 ml

1. Lyse cells (in 6 wells) with 500  $\mu$ l per well of Lysis buffer (+ protease inhibitors): 20 min at 4°C.
2. Clear by centrifugation (14000 rpm, 4°C 15 min Eppendorf centrifuge). Save an aliquot (30  $\mu$ l) for Western blotting.
3. Immunoprecipitate the kinase (e.g. with 10  $\mu$ l anti-flag affinity matrix beads, Sigma, for flag-tagged transfected kinase; or with appropriate antibody for endogenous kinase + Protein A-Sepharose or directly coupled to agarose): 2h at 4°C (rotating).
4. Wash the beads: 3x with 1 ml PBS (4°C), 1x with 1 ml Kinase buffer (4°C): pellet the beads by centrifugation (14000 rpm, 4°C, 45sec) and remove the supernatant.
5. Prepare Kinase buffer: add MnCl<sub>2</sub> to 10 mM (stock: 1 M) and <sup>32</sup>P- $\gamma$ -ATP (5  $\mu$ Ci per sample, usually 1/10 volume, i.e. 1  $\mu$ l of stock solution for one 10  $\mu$ l assay) and preincubate at 30°C for 10 min.
6. Add 1  $\mu$ g substrate: GST-I $\kappa$ B (1  $\mu$ l) or mutant substrate (as control) to the beads; add 10  $\mu$ l kinase buffer, mix gently and incubate at 30°C for 30 min (or longer).
7. Stop the reaction by addition of 4x SDS-sample buffer (4  $\mu$ l) and perform SDS-PAGE with the samples, followed by fixation of the gel (10% methanol, 10% HAc), drying and autoradiography.

For detection with PhastGel: use only 5  $\mu$ l beads, 5  $\mu$ l kinase buffer, 0.5  $\mu$ l substrate and 2  $\mu$ l 4x SDS-sample buffer: Run a 12.5% PhastGel with 4  $\mu$ l per sample