

	For final concentration of gel ( % T):					stack gel 5 ml
	Separating gel (10 ml)				Stack gel (10 ml)	
	7%	10%	12,5 %	15%	5%	
30% Acrylamide-bis solution 29:1 (A)	2.33	3.33	4.17	5	1.67	0.835
4x Separation buffer <b>1.5 M Tris/HCl pH 8.8</b>	2.5	2.5	2.5	2.5		0
4x Stacking buffer <b>0.5 M Tris/HCl pH 6.8 + phenol red</b>					2.5	1.25
aqua dest.	5	4	3.2	2.4	5.7	0
SDS (10 %)	0.1	0.1	0.1	0.1	0.1	2.85
TEMED	0.015	0.015	0.015	0.015	0.015	0.05
APS (10%)	0.03	0.03	0.03	0.03	0.03	0.0075

total ml	9.98	9.98	10.02	10.05	10.02	
Acryl in (A)	0.699	0.999	1.251	1.5	0.501	
Acryl %	7.01	10.02	12.49	14.93	5.00	

#### SDS-sample buffer

	final conc	6x	<b>6X Sample Buffer (with DTT)</b>	
Tris/HCl pH6.8	0.125	0.75	1 M Tris-Cl (pH 6.8)	2.4 ml
Glycerol	10%	60%	SDS	0.96 g
SDS	2%	12%	Glycerol	4.8 ml
DTT		739 mg	DTT	739 mg > ca. 600 m > final conc. 100mM
Bromophenol blue		4.8 mg	Bromophenol Blue	4.8 mg

Do not adjust the volume.

Reduction: either by 2%  $\beta$ -Mercaptoethanol  
(final conc., equals about 255 mM SH-groups, would mean about 125 mM DTT with 2 SH groups)

#### STRIPPING BUFFER (2% SDS, 62.5mM TRIS pH6.8, 100mM Beta-mercaptoethanol (bME))

Reagent	ml for 500ml	ml for 100 ml
20% SDS	50	10
1M TRIS pH6.8	31.25	6.25
H2O	418.75	83.75

#### \*\*ADD THE bME FRESH BEFORE

beta-Mercapto	3.55	0.71
---------------	------	------

Incubate the membrane for 15 - 30 min at 50°C submerged in stripping buffer.