

## **$\beta$ -Galactosidase Assay with CPRG**

- Lyse cells  
(recommended cell lysis buffer: 0.25M Tris/HCl pH 8.0, 0.25% (v/v) NP40, 2.5 mM EDTA)
- Pipet about 10  $\mu$ l of extract into a 96-well plate (appropriate neg. control / blank: 10  $\mu$ l of mock transfected cells, or non-transfected cells – as there is a slight endogenous  $\beta$ -Gal activity) – leave one well empty for blank (A-1)
- Add 100  $\mu$ l substrate solution to the wells (also to the blank-well)
- Incubate until red color develops (min to hours – depending on  $\beta$ -Gal activity, if you have low activity you can also incubate at 37°C)
- Optional: Stop with 50  $\mu$ l of Stop solution (only necessary if you want to time it exactly, e.g. by adding the substrate in a timed way and stopping the reaction in the same way)
- Measure with ELISA Reader at 570 nm (Filter #3)

### Lysis Buffer:

0.25M Tris/HCl pH 7.4 (or better 8.0)

0.25% (v/v) NP40

2.5 mM EDTA

CPRG-substrate solution: 1 mg/ml (= 1.65 mM)

in PBS + 10 mM KCl, + 1 mM MgCl<sub>2</sub>

alternative substrate buffer:

60 mM Na<sub>2</sub>HPO<sub>4</sub> pH 8.0, 1 mM MgCl<sub>2</sub>, 10 mM KCl, 50 mM Mercapto-ethanol

Stop solution: 0.5M Na<sub>2</sub>CO<sub>3</sub>